

Think. Inside. The. Box.

PURELAB® Prima 7/15/30*



Pre-treatment

LC140

Lifetime expectancy: 6 months

What is purified from the water?

Chlorine and particulates

Impact of not changing:

Chlorine, chloramine and organics remain, increasing the burden on downstream consumables which will need more frequent and expensive changes, and can damage RO



LC143

RO Membrane

Lifetime expectancy: 2-3 years

What is purified from the water?

Majority of organics and inorganics

Impact of not changing:

High levels of organics and inorganics remain in the water reducing the functional life of the DI pack



LC136M2

Vent Filter

Lifetime expectancy: 6 months

What protection does it provide?

Atmospheric filter, reducing CO₂ & airborne organics from entering the sealed system


Impact of not changing:

Increases in airborne contaminants and potential recontamination of purified water



*A combination of treatments and technologies produce the final water specification

For more information please contact your local ELGA water specialist:

 (877)-315 3542

 Elga.orders@veolia.com

 e-shop!

